



**The best  
modular cart  
solution for  
your application**



**MST 140: monitor stand**



The monitor stand is designed for the professional use of monitors and other equipment, such as laptops (with or without arms).

The monitor stand consists of the following parts; a frame with Ø 125mm. casters, two casters are fitted with a two-way brake and one mounting pole and side casings.

**Options:**  
A small shelf with handle (wxd: 460 x 435mm.), shelves without handle for front or back, a monitor arm with handle and VESA monitor adapter, bottle bracket, serum holder and scope holder.

**Specifications:**  
Height: 1400mm.  
Colour: RAL 9002 (off white).  
Max. load cart: 75kg. and 15kg. per shelf.

**SVR: open medical cart**



Open platform cart that is designed to make a compact set up for equipment, easily accessible from the front, sides and rear of the cart.

The basic SVR cart consists of the following parts; an open frame, including heavy-duty chrome plated Ø 125mm. casters, two casters are fitted with a two-way brake; a worktop and stainless steel poles, connecting the worktop to the base.

**Options:**  
Side casings, monitor rotating platform, closed bottom frame. Generic accessories e.g.: shelves, drawers, serum holder, camera holder, scope holder, lift- and extension arms, isolation transformer, medical grade monitor, motor drive system.

**Specifications:**  
Height: 1000mm., 1300mm., 1500mm. or 1700mm.  
Width: 460mm. or 540mm.  
Colour: RAL 9002 (off white).  
Max. load cart: 150kg. and 50kg. per shelf.

**MAC: secured medical cart**



Should you wish to secure the cart completely and keep the equipment safe the MAC can be equipped with double doors at the front. When opened the equipment is easily accessible.

The basic MAC cart consists of the following parts; closed bottom frame, including heavy-duty chrome plated Ø 125mm. casters, two casters are fitted with a two-way brake; a worktop; stainless steel poles, connecting the worktop to the base; side casings with handles and perforated steel rear door.

**Options:**  
Monitor rotating platform, generic accessories e.g.: shelves, drawers, serum holder, camera holder, scope holder, lift- and extension arm, isolation transformer, medical grade monitor, motor drive system.

**Specifications:**  
Height: 1300mm. or 1500mm.  
Width: 460mm. or 540mm.  
Colour: RAL 9002 (off white).  
Max. load cart: 150kg. and 50kg. per shelf.

**PS 130: multi use cart**



The PS 130 cart is an easy use cart that suites equipment up to 30kg. (TV, PC-monitor).

The basic cart consists of the following parts; an open frame with Ø 100mm. casters, two casters are fitted with a two-way brake; a worktop and stainless steel poles, connecting the worktop to the frame.

**Options:**  
Shelves for top frame, bottom frame and poles (adjustable in height), draw tray, mouse pad and drawers in two heights: 105mm. or 170mm.

**Specifications:**  
Height: adjustable from 950mm. to 1350mm. high.  
Width: 535mm. (nett. 460mm.)  
Colour: RAL 9002 (off white).  
Max. load cart: 80kg. and 50kg. per shelf.

**MDC: semi open medical cart**



The MDC is a platform cart with semi-open side walls. Equipment setup on this trolley is easily accessible from the front and rear. A perforated steel rear door shuts the cart at the rear, but still allows the equipment to be ventilated.

The basic MDC cart consists of the following parts; closed bottom frame, including heavy-duty chrome plated Ø 125mm. casters, two casters are fitted with a two-way brake; a worktop; stainless steel poles, connecting the worktop to the base; side casing with handles and a perforated steel rear door.

**Options:**  
Monitor rotating platform, generic accessories e.g.: shelves, drawers, serum holder, camera holder, scope holder, lift- and extension arms, isolation transformer, medical grade monitor, motor drive system.

**Specifications:**  
Height: 1300mm. or 1500mm.  
Width: 460mm. or 540mm.  
Colour: RAL 9002 (off white).  
Max. load cart: 150kg. and 50kg. per shelf.

**ECP 110: flexible endoscopy medical cart**



The ECP is a wide platform cart with semi open side walls that is designed to make a compact set up for your equipment that is easily accessible. The worktop is removable.

The basic ECP cart consists of a closed bottom frame, side casings with handles, heavy-duty chrome plated Ø 150mm. casters. A removable worktop and stainless steel poles.

**Options:**  
Set of doors at the rear, disposal container, fixed stand, adjustable stand, storage compartments and heavy-duty chrome plated Ø 200mm. casters. Two casters can be fitted with a two-way brake. To make it easier to move, two casters are fitted with fixed direction plates. Generic accessories e.g.: shelves, drawers, serum holder, camera holder, scope holder, lift- and extension arm, isolation transformer, medical grade monitor, motor drive system.

**Specifications:**  
Height: 1050mm.  
Width: 1190mm.  
Colour: RAL 9002 (off white).  
Max. load cart: 150kg. and 50kg. per shelf.



## Example of a basic cart



## Lift and extension arms



Jansen Medicars have been a specialist in the area of carts for medical applications since 1971. We design and produce platform carts that are used to make equipment for endoscopy and other medical applications mobile.

Their sturdy metal construction and a modular structure distinguish our products, enabling you to assemble a cart for any application. To meet specific wishes and requirements we have included additional products in our supply program, including lift / extension arm systems, medical grade monitors, isolation transformers and motor drive systems.

To assemble a cart you first select a type of cart. We have 6 basic carts: MST 140, PS 130, SVR, MDC, MAC and ECP (shown on previous page). Then you select the desired cart width and height and then the cart specific accessories. The cart can be completed with platforms, drawers, lift / extension arm(s), monitor(s), isolation transformer and a motor drive system.

For further options and details please consult our technical information sheet that can be send on request or downloaded from our website: [www.medicars.com](http://www.medicars.com).

## Shelves and drawers



## Generic accessories



## Motor drive system



## Electrical installation



## Example of a configured cart



[www.medicars.com](http://www.medicars.com)  
P.O. Box 49  
3600 AA Maarssen  
Netherlands

For over 35 years Jansen Medicars have designed, developed and produced furniture for medical equipment. Working closely with hospital specialists and suppliers of medical equipment allow us to produce ergonomically designed carts made from the best quality materials that meet the highest demands for medical and public use.

Since we design and produce the products ourselves, we are fully capable of producing custom work at attractive prices. Should you have specific wishes with respect to the products available, or if entirely new ideas are involved, we have the means to develop a product within a short period. Custom work is produced in close consultation with the client and as with all Jansen Medicars products, you receive a product that is durable and of the quality you are accustomed to.

Jansen Medicars also supply accessories and equipment to medical furniture such as isolation transformers, medical grade monitors, lift and extension monitor arms and motor drive systems.

Please feel free to visit our website: [www.medicars.com](http://www.medicars.com) or ask for our technical information sheet for more details.



**Dealer information:**

