


# NicoletOne™

MODULAR NEURODIAGNOSTIC SYSTEM



**VIASYS™**  
HEALTHCARE  
*Excellence For Life*

# NicoletOne™



*VIASYS Healthcare and Nicolet provide the highest quality products and service in the industry. NicoletOne is the perfect extension to any clinical practice.*

## SIMPLE, YET SOPHISTICATED

In today's healthcare environment, providing quality care requires a balance between the capabilities of diagnostic systems and their ease-of-use.

**NicoletOne™** is designed to achieve this balance by offering high-quality diagnostic information, in formats that are easy to use and customizable through an intuitive user interface **designed to shorten the learning curve.**

A key component of any neurodiagnostic application is the patient record and data management operation. Our **V-Link™** option **reduces errors and improves efficiency by integrating to your Hospital Information System (HIS).**

With a common patient/data management interface, application modules offering the level of stand-alone system functionality and a user interface **focused on minimizing training,** NicoletOne offers you a unique blend of **simplicity and sophistication.**



## MODULAR SCALABILITY

- One system can handle all your EEG, LTM, Sleep, EP, EMG and ICU testing needs.
- Add modalities as your demand requires.
- NicoletOne can improve your diagnostic capabilities and provide a backup to your primary neurodiagnostic systems.



## COMBINING PREMIER PRODUCTS

### EEG MODULE

The NicoletOne EEG module is a comprehensive, EEG system. Powerful, flexible and easy to use, the EEG module combines sophisticated signal analysis with innovative usability features making it the ideal choice for all of your EEG needs.



### SLEEP MODULE

Increase the functionality of your **EEG** or **LTM** module with the NicoletOne Sleep option. This option transforms your NicoletOne system into a fully-functioning Sleep machine, boasting a unique adaptable scoring method.

### VIKINGQUEST™ EMG MODULE

The VikingQuest EMG module brings the Nicolet Viking series' gold-standard in EMG and Nerve Conduction Studies to the NicoletOne system, from basic Motor Nerve studies all the way to Single Fiber EMG.

### LTM MODULE

The Long Term Monitoring (LTM) module offers unlimited, network tolerant recording with clear display of large montages. Unique compact digital amplifiers record 32 or 64 channels, or can be combined for up to 128 channels.

### VIKINGQUEST™ EVOKED POTENTIAL MODULE

The time-tested reliability of the VikingQuest system provides you with an Evoked Potential system for somatosensory, visual and auditory examinations you can depend on during diagnostic exams and surgery.

### MONITOR (ICU/nICU) MODULE

With touchscreen display, 16 or 32 channels and unlimited recording time, the NicoletOne Monitor is specifically designed with the challenges of the ICU and neonatal ICU environment in mind.

## PIONEERING CONNECTIVITY

- V-Link helps manage system connectivity, patient records, clinical data and workflow across all VIASYS NeuroCare products in your lab.
- All NicoletOne modules are fully integrated with the NicVue patient administrator; no need to learn a new interface.
- V-Link allows you to interface NicoletOne with your HIS.



## SUPERIOR FUNCTIONALITY

- Easy-to-use, high-quality amplifiers and stimulators result in the best data possible.
- Advanced software capabilities with user-friendly controls provide the answers you want – quickly.
- Leading interface technology available with NicoletOne is virtually impervious to interference from external sources.



## A CONTINUUM OF CONFIDENCE

Carrying on from its origins as Nicolet Biomedical, VIASYSHealthcare continues to **offer the leading technologies and the largest and best trained support organization** in the neurodiagnostics industry.

### COMPANY

With more than a century of combined experience, our trusted brands, Nicolet, Grason-Stadler and Toennies have provided **innovative technologies** for diagnostic audiology, neurodiagnostic, peripheral vascular and hand-held Doppler equipment.

### EXPERIENCED SUPPORT

Our more than 220 trained regional and in-house representatives provide **unmatched customer support**. We're committed to helping you get the most out of your VIASYS diagnostic and monitoring systems.

### HARDWARE

As with all our products, NicoletOne provides the perfect balance of **simplicity, ease of use, and flexibility**.

Quality amplifiers and precise measurement tools for each application provide complete reliability and accuracy.

### SUPPLIES

To complement the performance offered by all of our systems, we offer a full line of quality electrodes and supplies.



*VIASYS NeuroCare products continue to keep practices of all sizes connected with the best in neurodiagnostic technology.*



# NicoletOne™

**NicoletOne** is the modular neurodiagnostic system from VIASYS Healthcare. Whether you require a range of neurodiagnostic procedures or test protocols that dictate specific combinations of test modalities such as:

- EEG
- EMG
- NCV
- LTM
- EP,

the NicoletOne system is your answer.

## **SIMPLE, YET SOPHISTICATED**

NicoletOne's Windows-based user interface is simple to use, dramatically **speeding up the training process.**



## **COMBINING PREMIER PRODUCTS**


NicoletOne was conceived and designed to address the financial pressures facing healthcare professionals, such as:

- Managing patient volume
- Return on investment, and
- Future growth, without compromising the quality of care.

Each modular solution of the NicoletOne system is "**best in class,**" so you make no compromise with the quality of equipment you rely on.

**VIASYS Healthcare** – *Excellence for Life*





**VIASYS HEALTHCARE—EXCELLENCE FOR LIFE**

For more information about our audiology, neurodiagnostic, peripheral vascular and monitoring equipment, contact your VIASYS Healthcare professional. Offices in the U.S., Iceland, United Kingdom, Netherlands, Germany, and Japan with representative organizations worldwide.

P.O. BOX 44994 • MADISON WI 53744.4994 • PHONE: 608.273.5000 • TOLL FREE: 800.356.0007

GRASON-STADLER

NICOLET BIOMEDICAL

NICOLET VASCULAR

TOENNIES

[WWW.VIASYSHEALTHCARE.COM](http://WWW.VIASYSHEALTHCARE.COM)